



**THE ALOS-TTL
THROUGH THE LENS
VIDEO HEADLIGHT
SYSTEM**

**NOW WITH
VIVID IMAGE TECHNOLOGY**

QED Medical's ALOS Through the Lens Video Headlight System now brings line of sight viewing into crisp, HD focus with new Vivid Image Technology. View procedures through the eye's of the surgeon in the O.R., run live images to viewing rooms or record and edit for later use with incredible clarity.

Include Vivid Image Technology with your new QED Medical ALOS System or upgrade your existing QED System easily and affordably.

Integrated with the
Altair HI 2
Fiber Optic Headlight



Features:

- Brilliant, Clean, Even Field of Intensity
- Line of Sight Alignment with QED's TTL Design
- Variable Output including 720p and 1080p
- View in Real Time or Record
- Generous depth of field
- Operates with any QED Medical light source
- HDMI, BNC & S-Video Cables Included
- Available in NTSC or PAL
- 100v - 240v Power Supply

Full Systems Includes:*

- The classic OR level ALOS-TTL-XVI includes the**
- Altair HI 2 Fiber Optic Headlight w/1/4" CCD Camera
- and 15mm Lens
- HDVI Image Pack
- Nova 300-watt Xenon Light Source
- Camera Control Unit w/ Mounting Bracket
- System Stand
- HDMI, BNC/Composite & S-Video Cables

**Also available for use with Nova 200-watt Halogen and Beamer 10-watt LED Light Source and the C3 & XL Portable LED Headlight Systems.*

QED MEDICAL

SURGICAL LIGHTING PRODUCTS SINCE 1971

THE NAME BEHIND THE SYSTEMS YOU KNOW



US MADE SINCE 1971

WWW.QED-MEDICAL.COM
SALES@QED-MEDICAL.COM
750 ENTERPRISE DRIVE,
LEXINGTON, KY 40510
800-513-2256
859-231-0338 ph
859-231-0376 fax