

QED MEDICAL

US MADE SURGICAL LIGHTING SINCE 1971

the Beamer 10-watt LED Light Source

The Standard Beamer:

Small Footprint

Lightweight - only 10.5 Ounces

50,000 - hour LED

Runs Cool

3 - Port Turret Accommodates Popular Cable Ends

Earth Magnet Secures the Beamer to any Steel Surface

Provides Illumination for a Variety of Lighted Tools and Instruments

THE BEAMER can provide illumination for many different surgical tools, including:



Headlights



Lighted Retractors*



Endoscopes*

The Beamer - Portable:

Small, Clip-On Source

For Veterinary and other Mobile Applications

1.5 Hours on a Full Charge

3 Hour Charge Time

50,000 - hour LED

Single Port with several to choose from

Lightweight - only 13.5 Ounces

Provides Illumination for a Variety of Lighted Tools and Instruments

The Beamer:

QED-6001 Beamer 10-watt LED Light Source with Turret

QED-6025 Beamer Power Supply (included)

The Beamer - Portable

QED-6015 Beamer 10-watt LED Portable Clip On Single Port Source w/Battery


QED-6020 Battery Pack for QED-6015 (one included w/Light Source)

QED-6021 Charger for QED-6020 (one included w/ Light Source)



US MADE SINCE 1971

*These items not available through QED.



Beamer 10-watt LED Portable Clip On System

WWW.QED-MEDICAL.COM
SALES@QED-MEDICAL.COM
750 ENTERPRISE DRIVE
LEXINGTON, KY 40510
800-513-2256
859-231-0338 ph
859-231-0376 fax