

the C2 Portable Halogen Headlight System

Comfortable

Lightweight

Brilliant

Mobile

Adjustable

QED MEDICAL

US MADE SURGICAL LIGHTING SINCE 1971

Exceptional Intensity in a Portable Headlight

Ideal for ENT, Oral & Maxillofacial Surgery,
Plastic & Reconstructive Surgery and other Medical
and Dental applications

Untethered, Lightweight & Comfortable

100-Hour, 10-watt Halogen Lamp

2-Hour Battery Pack with Charger

Wall Transformer for Continuous Use

Convenient in-line On/Off Switch

Comfortable, Adjustable Headband

Clean, Even field of Intensity

Precision Spot Adjustment

Removeable, Sterilizable Positioning Bar

Full Motion Optic

The C2 System includes:

- QED-9500 C2, 10-watt Portable Halogen Headlight
on Standard Headband
- QED-9155 110v Wall Transformer*
- QED-9400 ECOPAC Rechargeable Battery Pack
- QED-9410 Battery Charger*

*Also Available for 220v - 240v Use

PORTABLE LAWSON LIGHT



***Also available with the Swivel Mount &
Monocular Optic ENT specific applications.***

Shown with
optional
Positioning
Bar

WWW.QED-MEDICAL.COM
SALES@QED-MEDICAL.COM
750 ENTERPRISE DRIVE
LEXINGTON, KY 40510
800-513-2256
859-231-0338 ph
859-231-0376 fax

US MADE SINCE 1971

