

# the Lawson Light Monocular Optic Headlight

Unique Design for ENT Applications

Swivel Positioning

Brilliant Intensity

Shadowless Field

# QED MEDICAL

US MADE SURGICAL LIGHTING SINCE 1971

Designed in cooperation with Dr. William Lawson  
of Mt. Sinai Medical Center, New York

Specifically created for the ENT Professional

Ball & Swivel Linkage allows the optic to be positioned  
over the dominate eye

A Tunnel drilled with precision through the optic  
provides unimpeded, co-axial, shadowless view into  
deep, restrictive areas of the body

Works with all QED Light Sources and others via  
the ACMI port

Single - Side Mount Cable:  
Lightweight & Flexible  
Reduces Neck Strain & Headlight Movement  
Easy & Affordable to Replace

Cables available in 7' & 10' Lengths with  
Black or Iridescent Blue Sheathing

Support & Flexibility built into every cable

Comfortable, Adjustable Headband

Collar & Gown Clips Included

Full Motion Optic

## ***The Lawson Light***

QED-7820 on Standard Headband w 7' Cable  
QED-7821 on Standard Headband w 10' Cable

***The Lawson Light provides exceptional intensity  
from any QED Medical Light Source and others.***

***Portable LED Lawson Light***



*Nova 300-watt Xenon*



*Nova 200-watt Halogen*



*Beamer 10-watt LED*

WWW.QED-MEDICAL.COM  
SALES@QED-MEDICAL.COM  
750 ENTERPRISE DRIVE  
LEXINGTON, KY 40510  
800-513-2256  
859-231-0338 ph  
859-231-0376 fax

US MADE SINCE 1971

