

# QED MEDICAL

US MADE SURGICAL LIGHTING SINCE 1971

## the Nova 3 200-Watt Halogen Light Source

Small Footprint

50-hour to 100-hour Lamp Life

3-Port Turret Accommodates Popular Cables

For Headlights & Lighted Instruments

White Light Filter Setting

Cooling Fan keeps the Nova 3 at Optimum Temperature

Easy Access Lamp Door

Handle for ease in movement

Lamp: 200 watt Halogen  
Lamp Life: 50 - 100 Hours  
Dimensions: 7.5" x 5.5" x 4.125"  
190mm x 139.7 mm x 104.7mm  
Weight: 6.75lbs / 3.06 kg  
Power: 100v - 240v 50/60 Hz  
Compliance: CAN/CSA C22.2 No. 601.1  
CAN/CSA C22.2 No. 601.1.2  
UL 60601-1  
IEC/EN 60601-1  
IEC/EN 60601-1-2

### **The Nova 3:**

QED-2050 200-watt Halogen Light Source w/Black Cover



Altair HI 2



Altair SA



Lawson Light



US MADE SINCE 1971



WWW.QED-MEDICAL.COM  
SALES@QED-MEDICAL.COM  
750 ENTERPRISE DRIVE  
LEXINGTON, KY 40510  
800-513-2256  
859-231-0338 ph  
859-231-0376 fax