

## THE NOVA 3

### 200-WATT HALOGEN LIGHT SOURCE

The Nova 3 packs brilliant halogen intensity into a small, feature laden footprint. Power any QED Medical headlight or other headlights & instruments via the three port turret.

Adjustable intensity, white light filter, thermostatically controlled fan and longer lamp life adjustments make the Nova 3 a powerful illumination option.



Optional System Stand



Altair HI 2

Lamp:	200 watt Halogen
Lamp Life:	50 - 100 Hours
Dimensions:	7.5" x 5.5" x 4.125" 190mm x 139.7 mm x 104.7mm
Weight:	6.75lbs / 3.06 kg
Power:	100v - 240v 50/60 Hz
Compliance:	CAN/CSA C22.2 No. 601.1 CAN/CSA C22.2 No. 601.1.2 UL 60601-1 IEC/EN 60601-1 IEC/EN 60601-1-2



Altair SA



Lawson Light  
Monocular Optic

The Nova 200-Watt Halogen Light Source powers each of these QED Medical Headlights, the ALOS-TTL Through the Lens Video Headlight System and, via the unique 3-Port Turret many other available headlights & instruments.

WWW.QED-MEDICAL.COM  
SALES@QED-MEDICAL.COM  
750 ENTERPRISE DRIVE,  
LEXINGTON, KY 40510  
800-513-2256  
859-231-0338 ph  
859-231-0376 fax

# QED MEDICAL

SURGICAL LIGHTING PRODUCTS SINCE 1971

THE NAME BEHIND THE SYSTEMS YOU KNOW



US MADE SINCE 1971