

QED MEDICAL

US MADE SURGICAL LIGHTING SINCE 1971

the Nova 300 300-watt Xenon Light Source

Brilliant, Variable Intensity

Color Correct Field

Simple to use

Immediate On/Off Function

1,000-hour Lamp

Dual Cooling Fans

5-Port Turret to Accommodate Major
Headlights

Known as a "Work Horse"
Easy to Change Lamps

Proven, Reliable Design

US Manufactured & Supported

Lamp: 300 watt Xenon
Color Temp: 6,000° Kelvin
Lamp Life: 1,000 Hours
Dimensions: 12.6" x 11.5" x 5.25"
320mm x 292mm x 133mm
Weight: 11 lbs / 5 kg
Power: 100v - 240v 50/60 Hz
Compliance: CAN/CSA C22.2 No. 601.1
CAN/CSA C22.2 No. 601.1.2
UL 60601-1
IEC/EN 60601-1
IEC/EN 60601-1-2

*QED has Private Labeled
these units for many partners.
Check the FDA Label on
the back to confirm your
unit is an original QED
Manufactured System.*



*Xenon remains the premier
source for consistent intensity
and color temperature*

 US MADE SINCE 1971



Let the Nova 300 Illuminate
these QED Medical or other Headlights



Altair HI 2



The Altair SA

The Nova 300:

QED-3010 Nova 300-watt Xenon Light Source

WWW.QED-MEDICAL.COM
SALES@QED-MEDICAL.COM
750 ENTERPRISE DRIVE
LEXINGTON, KY 40510
800-513-2256
859-231-0338 ph
859-231-0376 fax